# Prolaryn® Plus

## INJECTABLE IMPLANT

### VOCAL FOLD IMPLANT INSTRUCTIONS FOR USE

---

**Rx ONLY**

#### DESCRIPTION

Prolaryn® Plus injectable implant is a steam sterilized, latex-free, non-pyrogenic, semi-solid, cohesive, completely bio-degradable implant. The principal component of Prolaryn® Plus injectable implant is synthetic calcium hydroxylapatite, a biomaterial with over twenty years of use in orthopedics, neurosurgery, dentistry, otolaryngology and ophthalmology. Calcium hydroxylapatite is the primary mineral constituent of bone and teeth. The semi-solid nature of the implant is created by suspending calcium hydroxylapatite in a gel carrier that consists primarily of water (sterile water for injection USP) and glycerin (USP). The gel structure is formed by the addition of a small amount of sodium carboxymethylcellulose (USP). The gel is dissipated *in vivo* and replaced with soft tissue growth, while the calcium hydroxylapatite remains at the site of injection. The result is long-term yet non-permanent restoration and augmentation.

Prolaryn® Plus injectable implant is classified as a Class III Medical Device according to Annex IX of the MDD. Prolaryn® Plus injectable implant 1cc has a particle size range of 25-45 microns and can be injected with the two Rigid Injection needle(s) provided in this convenience pack. **For correct needle assembly please read the additional Information for Use, provided in this convenience pack, prior start of the treatment.**

#### INTENDED USE/INDICATIONS

Prolaryn® Plus injectable implant is indicated for vocal fold medialization and vocal fold insufficiency that may be improved by injection of a soft tissue bulking agent. The implant injection augments the size of the displaced or deformed vocal fold so that it may meet the opposing fold at the midline for improved phonation. Vocal fold insufficiency associated with serious aspiration difficulties may be an urgent indication.

#### CONTRAINDICATIONS

- Contraindicated in the presence of foreign bodies, acute inflammation, infection, inadequately controlled malignancy or rapidly advancing disease when these involve the laryngeal or upper respiratory tract.
- Contraindicated in bilateral laryngeal paralysis and vocal disorders of psychogenic or emotional origin.

#### WARNINGS

- Implant should not be injected into blood vessels. Injection into blood vessels may cause platelet aggregation, vascular occlusion, infarction, embolic phenomena or hemolysis.
- Airway obstruction following vocal fold injection can occur immediately or at anytime up to seven (7) days following injection. Airway obstruction results from aggressive vocal fold injection, over-injection, or laryngeal edema from trauma and manipulation of the larynx. Airway obstruction can often be prevented by intraoperative and postoperative steroid treatment and by minimizing laryngeal trauma and manipulation.

- Should not be injected into the airway. Confirm placement of needle tip visually before initiating implant injection.
- Should not be injected into organs or other structures that could be damaged by a space occupying implant.
- Some injectable implants have been associated with hardening of the tissues at an injection site, migration of particles from an injection site to other parts of the body and/or allergic or autoimmune reactions. Based on clinical usage, animal studies and supporting literature, this has not been observed nor is it expected with Prolaryn® Plus injectable implant.
- Safety and effectiveness during pregnancy or in lactating females has not been established.

## PRECAUTIONS

- Prolaryn® Plus injectable implant therapy should be delayed at least six months following the onset of the paralysis and/or until an adequate trial of voice rehabilitation has been made.
- Requires viable tissue for effectiveness. Scar tissue and significantly compromised tissue may not coapt appropriately. Patients who have had multiple surgeries may not be appropriate candidates for bulking injection treatment.
- Do not over-inject the vocal fold. Prolaryn® Plus injectable implant can be easily added in subsequent injections but cannot be easily removed. In extreme cases site rupture could occur.
- The Prolaryn® Plus injectable implant injection procedure and the associated instrumentation procedures have small but inherent risks of infection and/or bleeding like similar otolaryngology procedures. The patient may experience slight discomfort during and following the procedure. The usual precautions associated with otolaryngology procedures, specifically vocal fold injection, should be followed.
- Care should be taken with the injection of Prolaryn® Plus injectable implant, as with any surgical or implantation procedure, to avoid infection during the injection procedure. If infection occurs and cannot be corrected, it may become necessary to remove the implant.
- **Do not re-sterilize.** Prolaryn® Plus injectable implant is supplied sterile and non-pyrogenic in a sealed foil pouch and is intended for single patient, single treatment use only.
- The foil pouch should be carefully examined to verify that neither the pouch nor the syringe has been damaged during shipment. Do not use if the foil pouch is compromised or the syringe has been damaged. Do not use if the syringe end cap or syringe plunger is not in place. *There is a small amount of moisture normally present inside the foil pouch for sterilization purposes; this is not an indication of a defective product.*

## INCIDENT REPORTING

Any incident that directly or indirectly led, might have led or might lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, or a serious public health threat; and that has occurred in relation to the device should be reported under [AxCA-adverse.events@merz.com](mailto:AxCA-adverse.events@merz.com) or 1-877-336-4008 (Customer Service). In case of serious incident, please contact in addition [Ax-Safety@merz.de](mailto:Ax-Safety@merz.de). For Quality complaints please contact [AxCA-QA@merz.com](mailto:AxCA-QA@merz.com) or 1-877-336-4008 (Customer Service).

## POST MARKET SURVEILLANCE

The following adverse events have been identified during post-approval use of Prolaryn® Plus. Because they are reported voluntarily from a population of uncertain size, it is not always possible to reliably estimate their frequency or establish a causal connection to Prolaryn® Plus. These events have been chosen for inclusion due to a combination of their seriousness, frequency in reporting, or causal connection to Prolaryn® Plus: infection, over injection/under injection, ineffective treatment, tissue displacement, exposed material, superficial injection, edema, pain, stiffness at injection site, extraction of material, voice changes, dyspnea, airway obstruction, nodule, cough, hypersensitivity reactions including anaphylaxis, urticaria, erythema, facial edema, pruritus and hives.

## INDIVIDUALIZATION OF TREATMENT

Medical examination including medical history and diagnostic testing should be conducted to determine if the patient is an appropriate candidate for treatment with Prolaryn® Plus injectable implant. The outcome of treatment will vary between patients. In some patients additional treatments may be necessary to improve and/or maintain the level of resolution. If symptoms persist after treatment additional injections may be performed, but only after sufficient time has passed to evaluate the patient. The patient should not be re-injected sooner than seven days after the previous treatment.

## DIRECTIONS FOR USE

### PERCUTANEOUS VOCAL FOLD INJECTION

The following is required for the percutaneous vocal fold injection procedure:


- Prolaryn® Plus injectable implant filled syringe(s)
  - Appropriate size needle(s) {Provided separately}
  - Nasopharyngoscope
1. Prepare patient for nasopharyngoscopy and anesthetize using standard methods. Local anesthesia is not required but may be utilized at the injection site.
  2. Prepare the syringe(s), injection needle(s) and nasopharyngoscopic equipment before the surgical injection. A new injection needle may be used for each syringe, or the same injection needle may be connected to each new syringe.

Remove foil pouch from the carton. The pouch can be opened and the syringe dropped onto the sterile field when required. *There is a small amount of moisture normally present inside the foil pouch for sterilization purposes; this is not an indication of a defective product.*


Remove the Luer syringe cap from the distal end of the syringe prior to attaching the needle. The syringe can then be twisted onto the Luer lock fitting of the needle. **The needle must be tightened securely to the syringe and primed with Prolaryn® Plus injectable implant.** If excess Prolaryn® Plus injectable implant is on the surface of the Luer lock fittings, it will need to be wiped clean with sterile gauze. Slowly push the syringe plunger until Prolaryn® Plus injectable implant extrudes from the end of the injection needle. If leakage is noted at the Luer fitting, it may be necessary to remove the needle and clean the surfaces of the Luer fitting or, in extreme cases, replace both the syringe and the injection needle.

3. Place the nasopharyngoscope to precisely visualize the needle position and Prolaryn® Plus injectable implant volume during augmentation.

#### LOCATION AND DIRECTION OF NEEDLE PLACEMENT FOR TRANSCUTANEOUS INJECTION

[A] Through thyroid cartilage.

[B] Through cricothyroid membrane.

4. Extend the patient's neck if possible and identify by external landmark the cricoid and inferior border of the thyroid cartilage and thyroid notch. Because the superior surface of the vocal fold lies at approximately half the distance between the superior notch and the inferior border of the thyroid cartilage, injection is placed below this level but above the inferior thyroid cartilage margin. Transcartilaginous injection is used unless cartilage calcification prevents such injections, in which case needle placement is through the cricothyroid membrane.

**NOTE: Do not inject into a blood vessel.**

5. Prolaryn® Plus injectable implant should be injected lateral to the thyroarytenoid muscle. With needle location visually confirmed through the nasopharyngoscope, slowly push the plunger of the syringe to start the injection. After the initial injection, the patient should be asked to phonate and cough to disperse Prolaryn® Plus injectable implant throughout the vocal fold. Additional Prolaryn® Plus injectable implant is injected until the vocal folds touch during respiration at a position midway between the anterior commissure and the vocal processes.
6. Some tissue planes may be difficult to inject. If significant resistance is encountered when pushing the plunger, pull back the injection needle about one (1) to three (3) millimeters with the needle still in the vocal fold tissue. Push the plunger slowly again. If significant resistance is still encountered, it may be necessary to pull the needle entirely out of the injection site and try again in a new position. If significant resistance continues to persist, it may be necessary to try a different injection needle. If this is not successful, replace the syringe and injection needle.

#### ORAL VOCAL FOLD INJECTION

The following is required for the vocal fold injection procedure:

- Prolaryn® Plus injectable implant syringe(s)
- Appropriate size needle(s) {Provided separately}
- Nasopharyngoscope

1. Prepare patient for nasopharyngoscopy and anesthetize using standard methods. Local anesthesia is not required but may be utilized at the injection site.
2. Prepare the syringe(s), injection needle(s) and nasopharyngoscopic equipment before the surgical injection. A new injection needle may be used for each syringe or the same injection needle may be connected to each new syringe.

Remove foil pouch from the carton. The pouch can be opened and the syringe dropped onto the sterile field when required. *There is a small amount of moisture normally present inside the foil pouch for sterilization purposes; this is not an indication of a defective product.*

Remove the Luer syringe cap from the distal end of the syringe prior to attaching the needle. The syringe can then be twisted onto the Luer lock fitting of the needle. **The needle must be tightened securely to the syringe and primed with Prolaryn<sup>®</sup> Plus injectable implant.** If excess Prolaryn<sup>®</sup> Plus injectable implant is on the surface of the Luer lock fittings, it will need to be wiped clean with sterile gauze. Slowly push the syringe plunger until Prolaryn<sup>®</sup> Plus injectable implant extrudes from the end of the injection needle. If leakage is noted at the Luer fitting, it may be necessary to remove the needle and clean the surfaces of the Luer fitting or, in extreme cases, replace both the syringe and the injection needle.

3. Place the nasopharyngoscope to precisely visualize the needle position and Prolaryn<sup>®</sup> Plus injectable implant volume during augmentation.
4. Position injection needle for oral vocal fold injection. Prolaryn<sup>®</sup> Plus injectable implant should be injected lateral to the thyroarytenoid muscle.

**NOTE: Do not inject into a blood vessel.**

5. With needle location visually confirmed through the nasopharyngoscope, slowly push the plunger of the syringe to start the injection. After the initial injection, the patient should be asked to phonate and cough to disperse Prolaryn<sup>®</sup> Plus injectable implant throughout the vocal fold. Additional Prolaryn<sup>®</sup> Plus injectable implant is injected until the vocal folds touch during phonation at a position midway between the anterior commissure and the vocal processes.
6. Some tissue planes may be difficult to inject. If significant resistance is encountered when pushing the plunger, the injection needle should be pulled back about one (1) to three (3) millimeters (with the needle still in the vocal fold tissue) and push the plunger slowly again. If significant resistance is still encountered, it may be necessary to pull the needle entirely out of the injection site and try again in a new position. If significant resistance continues to persist, it may be necessary to try a different injection needle. If this is not successful, replace the syringe and injection needle.

#### **PATIENT COUNSELING INFORMATION**

- Provide an appropriate course of antibiotics as required.
- Instruct patient not to use the voice for three days. This minimizes any potential extrusion of the Prolaryn<sup>®</sup> Plus injectable implant through the injection site.

## HOW SUPPLIED

Prolaryn® Plus injectable implant is provided sterile (via steam) and non-pyrogenic in a syringe packaged in a foil pouch and boxed for convenient storage. Each syringe unit with needles convenience pack consists of one pre-filled syringe containing 1.0cc of Prolaryn® Plus injectable implant (degree of accuracy of syringe graduations is  $\pm 0.025\text{cc}$ ). Do not use if packaging and/or syringe are damaged or if the syringe end cap or syringe plunger is not intact. For correct needle assembly please read the additional Information for Use, provided in this convenience pack, prior start of the treatment.

**The contents of the syringe are intended for single patient, single treatment use only and cannot be re-sterilized. Re-use may compromise the functional properties of the device and/or lead to device failure. Re-use may also create a risk of contamination of the device and/or cause patient infection or cross-infection including but not limited to transmission of infectious disease(s) and blood transfer between patients. All which, in turn, may lead to patient injury, illness or death.**

## STORAGE

Packaged Prolaryn® Plus injectable implant should be stored at a controlled room temperature between 15°C and 32°C (59°F and 90°F). Do not use if the expiration date has been exceeded. The expiration date is printed on the product labels.

## DISPOSAL

Used and partially used syringes and injection needles could be biohazardous and should be handled and disposed of in accordance with facility medical practices and local, state or federal regulations.

## WARRANTY

Merz North America, Inc. warrants that reasonable care has been exercised in the design and manufacture of this product.

**THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ITS PARTICULAR PURPOSE.**

Handling and storage of this product, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Merz North America, Inc.'s control directly affect the product and the results obtained from its use. Merz North America, Inc.'s obligation under this warranty is limited to the replacement of this product and Merz North America, Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly, arising from the use of this product. Merz North America, Inc. neither assumes, nor authorizes any person to assume for Merz North America, Inc., any other or additional liability or responsibility in connection with this product.

Symbol	Title of symbol
	Catalogue number
	Batch code
	Unique device identifier
	Use-by date
	Do not re-sterilize
	Do not use if package is damaged
	Temperature limit
	Single Use Do not re-use
	Consult electronic instruction for use
	Sterilized using steam or dry heat
	Sterilized using ethylene oxide

 MANUFACTURED BY

Merz North America, Inc.  
4133 Courtney St., Suite 10  
Franksville, WI 53126 U.S.A.  
Telephone: 844.469.6379  
E-Mail: [mymerzsolutions@merz.com](mailto:mymerzsolutions@merz.com)

© 2024 Merz North America, Inc.  
Prolaryn is a registered trademark of Merz North America, Inc. The Merz logo is a trademark of Merz Pharma GmbH & Co. KGaA.

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© 2024, Merz North America, Inc.